

# Complete Drainage System

**SKATER™**



  
**ARGON**  
MEDICAL DEVICES

**The SKATER™ Drainage System** provides a complete portfolio of catheters and accessories for all-purpose, nephrostomy, biliary and centesis procedures.

### All-Purpose Drainage Catheters

For drainage of abscesses and fluid collections

- Choice Lock™ trocar stylet allows for single-handed abscess drainage placement
- Available in locking and non-locking configurations
- Expansive size offering of 6F - 16F
- Can be used safely with ethyl alcohol and doxycycline

### Kit Components Contain One:

- SKATER catheter
- Removable pigtail straightener
- Metal stiffening cannula
- Flexible stiffening cannula (Not available on 6F and 7F)
- Locking system with snap-off tab (Not available on non-locking pigtail catheters)

### Nephrostomy Drainage Catheters

For drainage of the urinary tract

- Designed for Seldinger technique
- Expansive size offering of 6F - 14F
- Can be used safely with ethyl alcohol and doxycycline

### Biliary Drainage Catheters

For drainage of the biliary tract

- Radiopaque marker band provides clear view of the most proximal drainage hole for accurate placement in the biliary tract
- Additional drainage holes for biliary duct to duodenum drainage
- Can be used safely with ethyl alcohol and doxycycline

