






















BLADES for Small-Bone Power Systems

Oscillating Saw Blades

MicroAire® Oscillating Saw (Models 1700, 1710, 1970, 2700, 2710, 5922, 6700, and 6710)

Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	
 CAUTION: Do not use the long 1200 Series sagittal blades in these models of oscillating saws, as they will overstress the oscillating mechanism.						
1700-001	10.2	11	0.5	0.3	36	
1700-002	10.2	12	0.5	0.3	36	
1700-003	13	16.5	0.5	0.3	36	
1700-007	10.2	11	0.5	0.3	36	
1700-008	10.2	12	0.5	0.3	36	
1700-009	10.2	7.7	0.5	0.3	36	
1700-010	10.2	9.3	0.5	0.3	36	
1700-012	10.2	10.2	0.5	0.3	36	
1700-013	9.1	7	0.5	0.3	36	 

MicroAire® and Hall® Oscillating Saw (Models 1740, 1972, 2740, 5972 and 6740)

Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	 CAUTION: These blades cannot be used in MicroAire® oscillating saw models 1700, 1710, 1970, 1975, 1976, 2700, 2710, 6700, and 6710.
ZO-060	5.8	12	0.6	0.4	36	
ZO-061	9.7	12	0.6	0.4	36	
ZO-062	5.8	18.5	0.6	0.4	36	
ZO-063	9.7	18.5	0.6	0.4	36	
ZO-064	5.8	28.5	0.6	0.4	36	
ZO-065	9.7	28.5	0.6	0.4	36	
ZO-066	18.5	28.5	0.6	0.4	36	
ZO-067	5.8	13.5	0.6	0.4	36	
ZO-068	9.7	13.5	0.6	0.4	36	
ZO-069	9.7	13.5	0.6	0.6	36	