


STERNUM SAW BLADES

Sternum Saw Blades

Linvatec®

Sternum Saw

Model 5059-05

Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	
ZR-031	9.3	53.6	1.0	0.64	18	

MicroAire®


Sternum Saw

Model 4800

Linvatec®

Sternum Saw


Model 5059-05


Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	
ZR-032	9.4	34	1.0	0.6	18	

MicroAire®

Sternum Saw

Model 7800

Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	
ZR-032M	9.4	34	1.0	0.6	18	

ZR-033M	9.4	34	1.0	0.6	16	
---------	-----	----	-----	-----	----	--

Linvatec®/ Hall®

Repeat Sternotomy Saw

Model 5057

Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	
ZO-079	53.3	30	0.6	0.4	18	