

Compression
Garment
Collection



CAROMED®
recover beautifully

CAROMED®

Head and Neck Garment Collection



Shaped Compression Vest

2-8002

2-8002OS

Black - Sizes: 34 - 44

Champagne - Sizes: 32 - 48

OS style only available in Champagne



2-8002



2-8002



2-8002OS

Features

- Provides comfortable compression to the chest
- Slightly contoured cup area
- Padded hook & eye closure behind front zipper for easy application
- 2" stabilization band placed low to avoid surgical site
- 2-2002 features front "Z" hook detachable straps
- 2-2002OS style features wide recloseable Velcro® shoulder straps

Ideal for

- Breast Reduction
- Breast Reconstruction
- Mastectomy
- Lumpectomy

Colors

Black

Champagne


Fabrics

Satin/Spandex

Post-Augmentation Bra

2-8124

Sizes: 26-54



Features

- Provides moderate compression to the chest
- Slightly contoured design for support and comfort
- Molded cups that shape and reduce pressure on the breasts
- Padded hook & eye beneath front zipper closure for easy application
- Break-away adjustable shoulder straps
- Seamless cup for smooth, uniform compression

Ideal for

- Breast Augmentation
- Breast Reduction
- Breast Reconstruction

Colors

Black

Champagne

Fabrics

Satin/Spandex

Unisex Multi-Flap Waistbinder

White:

3-8027 9" width

3-8037 12" width

Black:

3-8527 9" width

3-8537 12" width

3-8547 15" width

Sizes: XS - 3XL



3-8527-9



X-brace design for stability

Size	Waist	
	Inches	Centimeters
XS	22-24	55.9-61
S	26-28	66-71.1
M	30-36	76.2-91.4
L	38-42	96.5-106.7
XL	44-46	111.8-116.8
2XL	47-48	119.4-121.9
3XL	50-51	127-129.5

Features

- Provides maximum support and compression to the abdomen
- Multi-flap designs allows for greater contouring and adjustability
- Moisture wicking, cotton lining for comfort and smooth contact with skin
- Stacked horizontal Velcro® attachments for a greater range of adjustability and increased conformity
- X-brace design across back panel for stabilization
- Vertical side stays to minimize rollover
- Can be used with drainage tubing

Ideal for

- Abdominoplasty
- Liposuction of the abdomen and flanks
- Hernias

Colors

White

Black

Fabrics

Cotton/Nylon/Spandex