

7800

Cordless Sternum Saw

- 7800-002 Sternum Guard included
- 14.4-volt NiMH battery
- Three battery types
- Autoclavable or aseptic-transfer batteries
- Variable-speed trigger
- Optimal balance
- Simple twist-lock blade hub
- Standard or aggressive saw blade
- Adjustable for inferior or superior sternotomy approaches
- 12,000 CPM SV cycles per minute
- Handpiece and guard validated for use in automated washer/sanitizer

Accessories

ZR-032M	Sternum Saw Blade	p.18
7800-002	Sternum Guard	p.18
6640-710	NiMH Battery Pack	p.24
7505-710	Large NiMH Battery Pack	p.24
ZR-033M	Aggressive Saw Blade	p.18
7500-700	NiMH Battery Charger	p.24
7500-615	Aseptic Battery Housing	p.24
7500-620	Aseptic NiMH Battery Pack	p.24
7500-625	Charging Adapter	p.24
7500-630	Aseptic Transfer Shield	p.24



Part No.	Width	Usable Length	Cut Thickness	Material Thickness	TPI
ZR-032M	9.4 mm	34 mm	1.0 mm	0.6 mm	18
ZR-033M	9.4 mm	34 mm	1.0 mm	0.6 mm	16