

Cricothyrotomy Devices

Manujet III / Jet-Ventilation Catheters



The Manujet III with Jet-Ventilation Catheters is recommended for emergency rooms, crash carts, ambulances and operating rooms as it guarantees a quick and efficient oxygenation of a patient.

Emergency use

As a life saving manoeuvre in the “cannot intubate – cannot oxygenate” situation for oxygenation to avoid a severe desaturation of the patient.

For **pre-hospital use** if there is an obstruction of the upper airway. Trans Tracheal Jet-Ventilation (TTJV) is faster and simpler than surgical cricothyrotomy in emergencies with fewer complications.

Elective use

- ▶ Microlaryngoscopy
- ▶ Rigid bronchoscopy
- ▶ To assist a difficult fiberoptic intubation
- ▶ Predicted difficult extubation

Order information

Manujet III

With 4 m pressure hose, Jet-Ventilation Catheters acc. to Ravussin (13 G, 14 G, 16 G) and 100 cm connecting tube

REF 30-01-003

Manujet III

Ditto but with NIST plug and without pressure hose

REF 30-01-003NIST



Jet-Ventilation Catheters acc. to Ravussin

For single use, sterile

Size	Emergency Set with syringe	Box	Anaesthesia Set without syringe	Box
16 G	REF 30-02-918-1	1	REF 30-02-018-1	5
14 G	REF 30-02-914-1	1	REF 30-02-014-1	5
13 G	REF 30-02-913-1	1	REF 30-02-013-1	5

