



PREPARE TO BE IMPRESSED



impress[®]

Peripheral Angiographic Catheters

A revolutionary new design in
Peripheral Angiographic Catheters.



impress®

Peripheral Angiographic Catheters

Merit Medical introduces an impressive line of Peripheral Angiographic Catheters. By enhancing each aspect of the product from tip to hub, this catheter will not only meet your needs and expectations, it will impress you.

HUB

- Accordion design, insert-molded hub and strain relief is made of soft polyurethane material which minimizes kinking and provides steerability
- Ergonomically designed wings enhance handling and control
- The internal transition provides smooth guide wire advancement

