

ENDOCAVITY COVERS


NeoGuard, CIV-Flex and Eclipse covers are not made with natural rubber latex

- Viral barrier for protection against microbial migration, including viruses, bacteria, and bloodborne pathogens
- Recommended for transrectal and transvaginal procedures
- Sterile gel packets available separately



NeoGuard™ Covers

- Rolled material designed for ultrasound use
- Covers include colored elastic bands




CIV-Flex™ Covers

- Soft, durable material provides patient comfort
- Acoustically transparent, offering distortion-free imaging
- Slides easily over transducers
- Sterile cover kits include gel packet and colored elastic bands

Eclipse® Covers



- Pliable polyethylene material
- Pre-gelled with Aquasonic® 100 transmission gel
- Available for 2D and 3D/4D ultrasound transducers



REF #	COVER	PRODUCT DESCRIPTION	QTY.
610-1126	NeoGuard	Sterile 2 x 30cm (.8" x 11.8") cover	24
610-1127	NeoGuard	Non-sterile 2 x 30cm (.8" x 11.8") cover	24
610-843	NeoGuard	Sterile 2.6 x 30cm (1" x 11.8") cover	24
610-1038	NeoGuard	Sterile 2.6 x 30cm (1" x 11.8") cover	50
610-836	NeoGuard	Non-sterile 2.6 x 30cm (1" x 11.8") cover	24
610-1039	NeoGuard	Non-sterile 2.6 x 30cm (1" x 11.8") cover	50
676-117	NeoGuard	Non-sterile 2.6 x 30.5cm (1" x 12") cover for use with Endo-P II on G20 ultrasound system	24
610-844	NeoGuard	Sterile 4 x 30cm (1.6" x 11.8") cover	24
610-838	NeoGuard	Non-sterile 4 x 30cm (1.6" x 11.8") cover	24
610-792	Eclipse	Non-sterile 6cm tapered to 4.4 x 24.1cm (2.4" tapered to 1.8" x 9.5") polyethylene cover with gel	24
610-1102	Eclipse	Non-sterile 8.3cm tapered to 4.3 x 24.1cm (3.3" tapered to 1.7" x 9.5") polyethylene cover with gel	24
610-006	CIV-Flex	Sterile 11.9cm tapered to 3.8 x 61cm (4.7" tapered to 1.5" x 24") cover with gel	24
610-007	CIV-Flex	Non-sterile 11.9cm tapered to 3.8 x 61cm (4.7" tapered to 1.5" x 24") cover (no gel)	24

Latex Covers

- Loose-fitting, rolled cover
- Non-bulk covers conveniently packaged in individual perforated packets
- Cost-effective covers available in bulk-packaged, light-sensitive bags for extended shelf life



PROcovers™

- High-quality, custom-sized latex covers
- Covers include colored elastic bands

Tapered Latex Covers

- Provides form-fitting tip
- Offers easy application onto larger-handed transducers



REF #	COVER	PRODUCT DESCRIPTION	QTY.
610-773	Latex	Non-sterile 1.3 x 5.8cm (.5" x 2.3") cover for use with Series 500, 550	50
610-316	Latex	Non-sterile 1.5 x 7cm (.6" x 2.8") fingertip	100
610-774	Latex	Non-sterile 1.8 x 25.5cm (.7" x 10") cover with bands for use with Series 620	50
610-075	Latex	Sterile 2 x 20cm (.8" x 7.9") PROcovers	24
610-039	Latex	Non-sterile 2 x 20cm (.8" x 7.9") PROcovers	24
610-212	Latex	Sterile 2 x 30cm (.8" x 11.8") PROcovers	24
610-076	Latex	Sterile 2.6 x 20cm (1" x 7.9") PROcovers	24
610-040	Latex	Non-sterile 2.6 x 20cm (1" x 7.9") PROcovers	24
610-213	Latex	Sterile 2.6 x 30cm (1" x 11.8") PROcovers	24
610-101	Latex	Non-sterile 2.6 x 30cm (1" x 11.8") PROcovers	24
610-214	Latex	Sterile 3.5 x 20cm (1.4" x 7.9") cover with elastic bands	24
610-011	Latex	Non-sterile 3.5 x 20cm (1.4" x 7.9") cover with elastic bands	24
610-010	Latex	Non-sterile 3.5 x 20cm (1.4" x 7.9") cover with elastic bands	50
610-249	Latex	Non-sterile 3.5 x 20cm (1.4" x 7.9") cover - bulk	100
610-262	Latex	Non-sterile 3.5 x 20cm (1.4" x 7.9") cover - bulk	500
610-042	Latex	Sterile 4 x 30cm (1.6" x 11.8") PROcovers	24
610-041	Latex	Non-sterile 4 x 30cm (1.6" x 11.8") PROcovers	24
08990426	Latex	Sterile 3.2cm tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") cover with polyethylene cord cover, gel and elastic bands	10
610-985	Latex	Sterile 3.2cm tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") cover	24
610-795	Latex	Non-sterile 3.2cm tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") cover	24
08990418	Latex	Non-sterile 3.2cm tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") cover	100