

Enzymatic Detergents

Three Ways to Protect

Cleaning is the first and most important step in the sterilization process. These enzymatic detergents meet your needs in cleaning reusable medical instruments and devices. Our product offering is formulated to handle moderate to heavily soiled instruments and includes fragrance-free option for respiratory-sensitive individuals.



EmPower™ & EmPower™ Fragrance-Free Dual-Enzymatic Detergents

Formulated for moderate to heavily soiled instruments. Removes protein rich debris in both manual and automated reprocessing systems. Contains powerful cleaning agents that will not get cloudy at the use temperature.

Point-of-Use Products - EmPower Foam Enzymatic Spray & MetriSponge

The AAMI' Guideline on 'Safe Handling and Biological Decontamination of Reusable Medical Devices in Healthcare Facilities and in Non-Clinical Settings' recommends using products like EmPower Foam on instruments immediately after use to keep them moist in transport containers. MetriSponge is a contoured sponge saturated with extra-strength MetriZyme. Easily forms around endoscopes, the MetriSponge is a crucial part to removing soil which harbors bioburden before manual or automatic disinfection.

MetriZyme™ Dual-Enzymatic Detergent

Offers powerful cleaning, yet is gentle on instruments. MetriZyme is designed to perform effectively in the presence of high organic debris while offering broad cleaning action on a variety of protein soils. Surfactants provide additional cleaning on carbohydrates, lipids and protein soils. MetriZyme includes an enzyme stabilizer, a buffer to maintain pH, a preservative package, and solution dye to differentiate it from water and other solutions.

Ordering Information

Part No.	Description	Packaging
10-4102	EmPower 2-ounce	48 bottles/case
10-4100	EmPower 1-gallon	4 bottles/case
10-4150	EmPower 5-gallon	1 drum
10-4130	EmPower 30-gallon	1 drum
10-4400	EmPower Fragrance Free 1-gallon	4 bottles/case
NEW 10-4450	EmPower Fragrance Free 5-gallon	1 drum
10-4224	EmPower Foam - 24-oz. Spray	12 bottles/case
10-4025	MetriSponge	4 boxes/case
10-4005	MetriZyme - 1 quart	4 bottles/case
10-4000	MetriZyme - 1 gallon	2 bottles/case