

DRAIN

FEATURED PRODUCTS

- ReSolve® Biliary Locking Drainage Catheter
- ReSolve® Non-Locking Drainage Catheter
- ReSolve® Locking Drainage Catheter

**ReSolve® Biliary**  
Locking Drainage Catheter

**ReSolve®**  
Non-Locking Drainage Catheter

For all-purpose and biliary drainage procedures, ReSolve locking, non-locking, and biliary drainage catheters help to **LESSEN TRAUMA** and promote more accurate placement. The tip of the locking and biliary drainage catheters are located on the inside of the pigtail curve to prevent irritation of the surrounding tissue. The catheter shafts have depth-positioning markers to facilitate **ACCURATE** placement and to act as a reference for correct **POSITIONING**.

**ReSolve® Locking**  
Drainage Catheter

DRAIN