

VASCU-PICC® TLS

SMALL but POWERFUL

Approved for infusion therapy

From the smallest to the largest — and now everyone in-between — Medcomp® has a catheter size for any pediatric patient. VASCU-PICC® TLS is the latest addition to our PICC product line and features a dual lumen configuration, soft polyurethane body, and high flow rates in a small French size. Since smaller veins require greater precision, it only makes sense to offer you more choices in PICC sizes.

Especially with your smallest patients.

Comfort

Comfort is a special challenge with your pediatric patients. The smaller French size features a dual lumen which means fewer sticks, fewer tears.

Versatility

Since many pediatric patients require a combination of nutrition, fluids and/or medicine, the dual lumen configuration enables one catheter to serve multiple functions.

Perfect Fit

Medcomp® offers a full-line of pediatric PICC sizes to fit any pediatric patient. The new 2.6F fits nicely between existing 1.9F and 3F catheters, and makes catheter insertion easier.

2.6F Catheter Size 20cm & 50cm Lengths

Smaller catheter while maintaining dual access and reduced trauma in smaller vessels

Thermosensitive Material

Softens at body temperature

Easy-to-read Depth marks

Multiple Function Luers

Eliminates separate access lines and decreases vessel trauma

Dual Lumen Design

- Combined access to one site reduces chance of infection, providing better patient outcomes and increased parent satisfaction
- Offers the ability to infuse incompatible drugs

Polyurethane Material

- Provides strength, versatility, and high flow rates in a smaller French size
- Higher flow rates and less chance of occlusion
- Allows greater control during insertion

Radiopaque Tip identifiable under X-Ray

