


Infu-Stat™

Disposable Pressure Infuser



Infu-Stat™ — A pressure infuser bag designed for single patient use to prevent cross contamination and eliminate costly cleaning and repair. The Infu-Stat™ gauge incorporates a pressure release at 325mm/hg and includes large easy-to-read numbers.

**20 Years of Medical
Device Sales**

 **Mason
Tayer**
Medical Products Corp.

Disposable Pressure Infuser

Infu-Stat™

Economical, single patient use to prevent cross contamination and eliminate costly cleaning and repair.

- ◆ **Loads easily.**
- ◆ **One hand operation** for inflation or deflation.
- ◆ **Large bore stopcock** for rapid inflation and deflation.
- ◆ **Large numbers & color-coded graphics** for easy gauge pressure reading.
- ◆ **Gauge** maintains reliable pressure.
- ◆ **Safety valve** eliminates the possibility of over inflation and releases air at 325mm/hg.
- ◆ **3 different volume bags** (500ml, 1000ml, 3000ml) color-coded for easy identification.
- ◆ **Latex free.**



Item#	Size	Packaging
MTM305	500ml	5/box - 25/cs
MTM310	1000ml	5/box - 25/cs
MTM330E	3000ml	2/box - 10/cs

 **OVERNIGHT SHIPPING**

 **INTERNATIONAL SHIPPING**

 **MINIMAL PACKAGING**

We gladly accept:



For additional information please visit
www.MasonTayler.com



Toll Free US and Canada: 1-888-270-7494
International: 001-716-883-8097
Fax: 716-883-8099
Email: info@masontayler.com
1083 Delaware Ave., Buffalo, NY 14209 USA

Made in USA