

SAFEGUARD RADIAL™

why safeguard radial?

The Safeguard Radial Compression Device is designed to help you achieve **PATENT HEMOSTASIS** at the conclusion of your radial procedure. The device is **COMFORTABLE** for your patients, **CONVENIENT** for hospital staff, and **EASY** to use.

> patent hemostasis


- Focuses compression only on the radial artery where you want it
- Puncture site clearly visible for visual assessment
- Stable positioning of bulb for optimal compression, not occlusion

> comfort

- Non-woven textile band for patient comfort
- No rigid plastic edges
- Low profile design

> convenient and easy

- One size fits all, easily adjusted to any wrist size
- Uses standard Luer lock syringe, no need for special syringe
- Easy to monitor distal and proximal pulses to ensure patent hemostasis is achieved.
- Easy inflation and deflation



SAFEGUARD RADIAL™
Compression Device

MERIT MEDICAL RADIAL™

THINK RADIAL. DISCOVER MERIT.™

ORDER INFORMATION

Catalog Number:

Description:

Length

SGR26 (Former catalog number 0684-00-0572)

SAFEGUARD Radial Compression Device

26 cm

Quantity: 10 units per box, packaged individually
Contact Merit Medical for more information.