



The ergonomic Solo™ Driver is the smallest wire driver with the widest range of wire sizes, part of the new 5000 Series small-bone power system developed with surgeons, for surgeons by MicroAire.

FEATURES AND BENEFITS

- Single-hand operation
- Fingertip control (forward, reverse, oscillate)
- Optimal visibility for hand and foot procedures
- Fully cannulated
- Light Weight and compact
- Ergonomic - resulting in less hand fatigue for surgeons
- No-fail wire guidance system
- Superior control for small delicate applications
- Ambidextrous

Every virtuoso requires the right instrument.



5995U
Solo™ Driver
(Unitized) Wire Driver

The ergonomic Solo™ Driver is the smallest wire driver with the widest range of wire sizes, part of the 5000 Series small-bone power system developed with surgeons, for surgeons (for use with the new 5025 electric power console).

- Single-hand operation
- Fingertip control (forward, reverse, oscillate)
- Optimal visibility for hand and foot procedures
- Fully cannulated
- Wide range of wires (0.7mm to 2.0mm)
- Smaller than other wire-drivers
- Designed for ambidextrous use
- Ergonomic, comfortable and powerful

Accessories		
5025	Electric Instrument Control Console	p.25
5006-5000	Electric Instrument Cable for 5025 Console	p. 25
View K-Wires & Pins in the Power Disposables Catalog at www.microaire.com		

