



SURGICAL COVERS

NeoGuard™, CIV-Flex™ and Polyethylene covers are not made with natural rubber latex

- Viral barrier for protection against microbial migration, including viruses, bacteria, and bloodborne pathogens
- Extends field of sterile protection against cross-contamination
- Accordion- or telescopically-folded for easy application; Flat-folded cover also available

Surgi™ Covers


REF #	FOLD	PRODUCT DESCRIPTION	QTY.
610-367	Accordion	Sterile 15.2cm tapered to 7.6 x 244cm (6" tapered to 3" x 96") polyethylene cover with 14 x 30cm (5.5" x 12") CIV-Flex cover	12
610-023	Accordion	Sterile 15.2cm tapered to 7.6 x 244cm (6" tapered to 3" x 96") polyethylene cover with 14 x 30cm (5.5" x 12") CIV-Flex cover (no gel)	12
610-217	Flat	Sterile 15.2cm tapered to 7.6 x 244cm (6" tapered to 3" x 96") polyethylene cover with 14 x 30cm (5.5" x 12") CIV-Flex cover	24
610-664	Accordion	Sterile 15.2cm tapered to 7.6 x 244cm (6" tapered to 3" x 96") polyethylene cover with nose cone (small) for use with Philips (ATL Ultrasound) Access A, B, C, 725B, WAAA transducers	10
610-665	Accordion	Sterile 15.2cm tapered to 7.6 x 244cm (6" tapered to 3" x 96") polyethylene cover with nose cone (large) for use with Philips (ATL Ultrasound) Access A, B, C, 725B, WAAA transducers	10



NeoGuard™ Surgi-Tip™ Covers

- Provides snug fit to small-footprint transducers
- Cover includes sterile gel packet, probe clips and colored elastic bands

REF #	FOLD	PRODUCT DESCRIPTION	QTY.
610-1072	Accordion	Sterile 15.2 x 244cm (6" x 96") polyethylene cover with attached 3cm (1.2") NeoGuard tip	12
610-833	Telescopically	Sterile 15.2 x 244cm (6" x 96") polyethylene cover with attached 3cm (1.2") NeoGuard tip	12
610-834	Telescopically	Sterile 20.3 x 244cm (8" x 96") polyethylene cover with attached 6cm (2.4") NeoGuard tip	12



NeoGuard™ Surgi-Boot™ Cover


- For use with Esaote Omni, GE Healthcare i12L, Philips (ATL Ultrasound) CL10-5, CL15-7, L15-7io, Siemens VF13-5SP, SonoSite HST, SLA/13-6, Supersonic SHL20-6 transducers
- Form-fitting tip provides snug fit to small footprint transducers

REF #	FOLD	PRODUCT DESCRIPTION	QTY.
610-797	Telescopically	Sterile 15.2 x 244cm (6" x 96") polyethylene cover with attached form-fitting NeoGuard tip	12



T-Shaped Surgi™ Cover

REF #	FOLD	PRODUCT DESCRIPTION	QTY.
645-058	Accordion	Sterile 15.3 x 244cm (6" x 96") polyethylene cover with attached form-fitting latex tip	12



ACUSON AcuNav™ Cover

- For use with ACUSON AcuNav ultrasound catheter family connector

REF #	FOLD	PRODUCT DESCRIPTION	QTY.
612-113	Accordion	Sterile 15.2cm tapered to 7.6 x 183cm (6" tapered to 3" x 72") polyethylene cover	12



Biosense Soundstar® eco Catheter Connector Cover


REF #	FOLD	PRODUCT DESCRIPTION	QTY.
608-001	Accordion	Sterile 20.3cm tapered to 10.8 x 183cm (8" tapered to 4.3" x 72") polyethylene cover	12



Laparoscopic Cover

- Form-fitting tipped transducer cover fits standard 10mm laparoscopic transducers (requires 12mm trocar)
- Cover is not made with natural rubber latex

REF #	FOLD	PRODUCT DESCRIPTION	QTY.
610-941	Telescopically	Sterile 15.2 x 244cm (6" x 96") polyethylene cover with attached form-fitting tip, clips	12



Cord Cover & System Drapes

- Probe clips and colored elastic bands included with cord covers
- Adhesive tape included with system drapes


REF #	PRODUCT DESCRIPTION	QTY.
610-024	Sterile 7.6 x 244cm (3" x 96") polyethylene cord cover	24
610-1241	Sterile 27.9 x 32.4cm (11" x 12.75") polyethylene system drape, tape strips (3)	24
610-1037	Sterile 50.8 x 71.1cm (20" x 28") polyethylene system drape, tape strips (4)	24
610-025	Sterile 76 x 76cm (30" x 30") polyethylene system drape	24
DR-000537	Sterile polyethylene instrument drape for use with NeuWave Certus™ 140 PDM	24



Viewing Wand Cover

- For use with ISG Technologies Viewing Wand

REF #	PRODUCT DESCRIPTION	QTY.
650-023	Sterile 30.5 x 117cm (12" x 46.1") polyethylene wand cover	30



TEE Covers

- Covers are not made with natural rubber latex
- Tapered for tighter fit on TEE transducers

REF #	PRODUCT DESCRIPTION	QTY.
610-933	Sterile 4.6cm tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex™ cover on applicator, gel packet	12
610-840	Non-sterile 4.6cm tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe, bite block	12
610-860	Non-sterile 4.6cm tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator	12
610-1013	Sterile 3.8cm tapered to 1.7 x 91.4cm (1.5" tapered to .7" x 36") CIV-Flex cover on applicator, gel packet	12
610-1014	Non-sterile 3.8cm tapered to 1.7 x 91.4cm (1.5" tapered to .7" x 36") CIV-Flex cover on applicator, gel-filled syringe, bite block	12
610-1229	Non-sterile 3.8cm tapered to 1.7 x 91.4cm (1.5" tapered to .7" x 36") CIV-Flex cover on applicator	12



C-Arm Drapes

- C-Arm drapes are quick to set up, efficient in use, and easy to discard
- Drapes are individually pouched for single-use
- Allow for equipment and patient protection during procedures
- Drapes made from durable polyethylene to provide a proper fit for any surgical procedure



REF #	PRODUCT DESCRIPTION	QTY.
610-5429	Sterile 9" / 12" I.I. C-Arm kit includes C-Arm drape with 5 mounting clips, I.I. cover, X-Ray tube cover, footswitch cover fits OEC 9600/9800, Siemens ARCADIS Avantec, Siemens SIREMOBIL Compact L, Siemens ARCADIS Varic, Philips BV series, Ziehm	20
610-5411	Sterile 9" I.I. / X-Ray tube cover, 30" x 32" contoured bandbag (<i>individually pouched</i>)	60
610-5412	Sterile 12" I.I. / X-Ray Tube cover, 36" x 36" contoured bandbag (<i>individually pouched</i>)	30
610-5414	Sterile 9" I.I. universal halfbag cover, 36" x 80" designed to fit I.I. and midway down C-Arm	20
610-5434	Sterile Mobile C-Arm drape, clear poly 42" x 74" drape with poly straps (<i>individually pouched</i>)	20
610-5402	Sterile Mini C-Arm kit with MiniView C-Arm drape, keyboard cover, footswitch cover fits OEC Mini 6800/6600 (<i>individually pouched</i>)	20
610-5403	Sterile Mini C-Arm drape - single piece drape to cover I.I., X-Ray tube and C-Arm fits OEC Mini 6800/6600 (<i>individually pouched</i>)	25
610-5423	Sterile Mini C-Arm drape - single piece drape to cover I.I., X-Ray tube and C-Arm fits Fluoroscanner /Hologic (<i>individually pouched</i>)	20

Urology Disposables



REF #	PRODUCT DESCRIPTION	QTY.
610-5311	Sterile Uro-Catcher drain bag fits all hoop style frames	10
610-5313	Sterile Uro-Catcher drain bag fits all hoop style frames	20
610-5416	Sterile Urology Disposables Pack with 5" x 10" hand control cover, 24" x 20" footswitch cover, 20" x 36" footswitch cover fits OEC UroView 2600/2800	20
610-5442	Sterile Uro-Catcher drain bag fits Siemens Uroskop System and Dornier	20

Miscellaneous Disposables

Bands are not made with natural rubber latex

REF #	PRODUCT DESCRIPTION	QTY.
610-655C	Sterile 36" x 30" band bags, clear polyethylene, elastic band (<i>individually pouched</i>)	100
610-657C	Sterile 36" x 36" band bags, clear polyethylene, elastic band (<i>individually pouched</i>)	100
610-668C	Sterile 42" x 36" band bags, clear polyethylene, elastic band (<i>individually pouched</i>)	100
610-704C	Sterile 54" x 54" slush drape, polyurethane (<i>individually pouched</i>)	25
610-705C	Sterile 24" x 36" X-Ray cassette cover (<i>individually pouched</i>)	25
610-733C	Sterile 48" (<i>expanded width</i>) x 24" shower cap style bag, elastic opening	50

