




Help Reduce Medication Errors

MERIT MEDALLION® SERIES SYRINGES



*Merit Medallion Syringes
help take the guess work out of
medication labeling.*



Color coded, pre-labeled syringes can help reduce medication errors



The Medallion Syringes are available in :

7 Colors: Light Green, Purple, Yellow, Red, Dark Green, Light Blue and White

Multiple Sizes:
1,2,3,6,10,20,30 & 60ml

Customized Printing:
All medication names can be printed directly on the syringe barrel in white or black ink.

Luer Options:
Fixed Male Luer or Slip Luer



Polycarbonate syringe barrels create a “smooth as glass” feel and help increase the visualization of fluid and air bubbles.



Custom Medallion syringes can include a unique frosted field for doctors or clinicians to write additional pertinent information.



Medallion® Syringes have clear polycarbonate barrels, optional colored plungers, customized medication printing and fixed male or slip luers.

Select one option from each category and identify your custom medication to be printed on the barrel.

Your Merit Representative will provide you with the corresponding catalog number.

Custom Medallion® Syringes are packaged 25 units per box, 4 boxes per case.

NOTE: Custom printed Medallion® Syringes are available upon special request. Please allow up to six weeks for delivery.



The combination of color coding and custom printed syringes help prevent potentially dangerous medication mix-ups.



The Medallion's high quality printing will not smudge or peel away from the syringe barrel.



JCAHO
PATIENT SAFETY GOALS
"Labeling occurs when any medication or solution is transferred from the original packaging to another container."

*Source: Joint Commission
 2009 National Patient
 Safety Goals Elements of
 Performance for NPSG.03.04.01*



Rise above the rest with color coded Medallion Syringes.