

# ENDOCAVITY NEEDLE GUIDES

CIVCO's endocavity needle guidance systems are designed to correspond with on-screen software guidelines, directing instruments for accurate placement procedures. Manufactured for a variety of transducers, the needle guides are designed to provide a secure fit to transducers. CIVCO needle guides are available in both disposable and reusable configurations.

## Endocavity Needle Guidance Solutions

CIVCO's complete line of endocavity needle guidance solutions include transrectal and transvaginal guides for ultrasound transducers from leading ultrasound system manufacturers and are recommended for procedures such as: tissue biopsies, fluid aspirations, catheterization and marker placement.



**Disposable Needle Guides** reduce risks associated with cross-contamination. According to the Centers for Disease Control and Prevention, “hospital-acquired infections account for an estimated 1.7 million infections, 99,000 deaths and \$30 billion in excess healthcare costs annually.” By using a disposable guide and a transducer cover, you can help fight infection in your facility.

The guides are conveniently packaged in sterile procedure kits containing needle guides, transducer cover, gel packet and colored elastic bands. Most guides accept 16-18GA instruments.



**Reusable Needle Guides** are packaged non-sterile. The reusable guides require proper sterilization and disinfection after each use. Most guides accept 16-18GA instruments.



**Disposable Quick-Release Needle Guides** are designed for single-use and provide physicians with a tool for performing procedures such as: deep perirectal abscess drainage, core tissue biopsy and larger-gauge device and catheter placement. The quick-release device offers three sizes in a single convenient procedure kit. Gauges include: 10-11GA, 14GA and 16-18GA.

For cleaning brushes and other endocavity accessories, visit [civco.com](http://civco.com)